

Product Datasheet

SARS-CoV-2 NSP9 polyclonal antibody (orb1718881)

Description	SARS-CoV-2 NSP9 polyclonal antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of coronavirus NSP9(YP_009742616.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	15-20KDa
Uniprot ID	P0DTD1
Dilution Range	WB, 1:500 - 1:2000
Expiration Date	12 months from date of receipt.