

Product Datasheet

ZCCHC7 polyclonal antibody (orb1719219)

Description	ZCCHC7 polyclonal antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ZCCHC7(NP_115602.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	63kDa
Uniprot ID	Q8N3Z6
Dilution Range	WB, 1:500 - 1:2000
Expiration Date	12 months from date of receipt.