

Product Datasheet

Neuropilin-1 (NRP1) polyclonal antibody (orb1719952)

Description	Neuropilin-1 (NRP1) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic peptide of human Neuropilin-1 (NRP1)(NP_003864.5).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	120KDa
Uniprot ID	O14786
Dilution Range	WB,1:500 - 1:3000 IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.