

## Product Datasheet

### SIRT4 polyclonal antibody (orb1720379)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | SIRT4 polyclonal antibody                                                                |
| <b>Species/Host</b>        | Rabbit                                                                                   |
| <b>Reactivity</b>          | Human                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | IF                                                                                       |
| <b>Immunogen</b>           | Recombinant fusion protein of human SIRT4(NP_036372.1).                                  |
| <b>Preservatives</b>       | 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2                              |
| <b>Concentration</b>       | 1mg/ml                                                                                   |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                               |
| <b>Uniprot ID</b>          | <b>Q9Y6E7</b>                                                                            |
| <b>Dilution Range</b>      | IF/ICC, 1:50 - 1:200                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |