

Product Datasheet

VPS26B polyclonal antibody (orb1720682)

Description	VPS26B polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein of human VPS26B(NP_443107.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q4G0F5
Dilution Range	IF/ICC, 1:50 - 1:200
Expiration Date	12 months from date of receipt.