

Product Datasheet

ZSWIM1 polyclonal antibody (orb1720815)

Description	ZSWIM1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ZSWIM1(NP_542170.3).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	55kDa
Uniprot ID	Q9BR11
Dilution Range	WB, 1:500 - 1:2000
Expiration Date	12 months from date of receipt.