

Product Datasheet

ITGB3BP polyclonal antibody (orb1720866)

Description	ITGB3BP polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human ITGB3BP(NP_055103.3).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q13352
Dilution Range	IHC,1:50 - 1:200
Expiration Date	12 months from date of receipt.