

Product Datasheet

ATP6V1A polyclonal antibody (orb1720951)

Description	ATP6V1A polyclonal antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human ATP6V1A(NP_001681.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	74kDa
Uniprot ID	P38606
Dilution Range	WB, 1:500 - 1:2000
Expiration Date	12 months from date of receipt.