

Product Datasheet

CLEC4C polyclonal antibody (orb1721441)

Description	CLEC4C polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CLEC4C(NP_569708.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	38kDa
Uniprot ID	Q8WTT0
Dilution Range	WB, 1:500 - 1:2000
Expiration Date	12 months from date of receipt.