

Product Datasheet

TNFR2/TNFRSF1B polyclonal antibody (orb1721489)

Description	TNFR2/TNFRSF1B polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TNFR2/TNFR2/TNFRSF1B(NP_001057.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	60kDa
Uniprot ID	P20333
Dilution Range	WB, 1:500 - 1:2000
Expiration Date	12 months from date of receipt.