

Product Datasheet

MDH2 polyclonal antibody (orb1721503)

Description	MDH2 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MDH2(NP_005909.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	35kDa
Uniprot ID	P40926
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.