

Product Datasheet

PACS1 polyclonal antibody (orb1721901)

Description	PACS1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PACS1(NP_060496.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	105kDa
Uniprot ID	Q6VY07
Dilution Range	WB, 1:1000 - 1:3000
Expiration Date	12 months from date of receipt.