

Product Datasheet

KDM4A polyclonal antibody (orb1721954)

Description	KDM4A polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human KDM4A(NP_055478.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Refer to figures
Uniprot ID	O75164
Dilution Range	IHC,1:50 - 1:200
Expiration Date	12 months from date of receipt.