

Product Datasheet

RNF4 polyclonal antibody (orb1722142)

Description	RNF4 polyclonal antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human RNF4(NP_001171938.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	21kDa
Uniprot ID	P78317
Dilution Range	WB, 1:500 - 1:1000
Expiration Date	12 months from date of receipt.