

Product Datasheet

APEX1/APE1 polyclonal antibody (orb1722234)

Description	APEX1/APE1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human APEX1/APE1(NP_542380.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	39kDa
Uniprot ID	P27695
Dilution Range	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:100
Expiration Date	12 months from date of receipt.