

Product Datasheet

Niemann-Pick type C1 Like-1 polyclonal antibody (orb1722523)

Description	Niemann-Pick type C1 Like-1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Niemann-Pick type C1 Like-1(NP_001095118.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	148kDa
Uniprot ID	Q9UHC9
Dilution Range	WB,1:1000 - 1:4000
Expiration Date	12 months from date of receipt.