

Product Datasheet

Lactic acid-Histone H3-K14 polyclonal antibody (orb1722737)

Description	Lactic acid-Histone H3-K14 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic lactylated peptide around K14 of human Histone H3(NP_003520.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	17KDa
Uniprot ID	P68431
Dilution Range	WB, 1:500 - 1:2000
Expiration Date	12 months from date of receipt.