

Product Datasheet

Pan DiMethyl-lysine polyclonal antibody (orb1722739)

Description	Pan DiMethyl-lysine polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	A synthetic peptide corresponding to a sequence containing dimethylatedK.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	15-60KDa
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:100 IF/ICC,1:50 - 1:100 IP,1:50 - 1:100
Expiration Date	12 months from date of receipt.