

Product Datasheet

MonoMethyl-Histone H3-K36 polyclonal antibody (orb1722740)

Description	MonoMethyl-Histone H3-K36 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, IP, WB
Immunogen	A synthetic monomethylated peptide around K36 of human histone H3(NP_003520.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	17KDa
Uniprot ID	P68431
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:200 IF/ICC,1:50 - 1:200 IP,1:50 - 1:200 ChIP,1:20 - 1:100 ChIP-seq,1:20 - 1:100
Expiration Date	12 months from date of receipt.