

Product Datasheet

DiMethyl-DNMT3A-K44 polyclonal antibody (orb1722750)

Description	DiMethyl-DNMT3A-K44 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	A synthetic dimethylated peptide around K44 of mouse DNMT3A(NP_031898.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O88508
Dilution Range	IF/ICC, 1:50 - 1:200
Expiration Date	12 months from date of receipt.