

Product Datasheet

Annexin A1 monoclonal antibody (orb1722757)

Description	Annexin A1 monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human Annexin A1 beta protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Uniprot ID	P04083
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.