

Product Datasheet

Rad9 monoclonal antibody (orb1722762)

Description	Rad9 monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Purified recombinant human RAD9A protein fragments expressed in E.coli
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 43 kDa; Observed MW: 55 kDa
Uniprot ID	Q99638
Dilution Range	WB: 1/500-1/1000 IP: 1/20
Expiration Date	12 months from date of receipt.