

Product Datasheet

Hydroxyl-Histone H2A (Tyr39) monoclonal antibody (orb1722780)

Description	Hydroxyl-Histone H2A (Tyr39) monoclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A synthesized peptide derived from human Histone H2A (hydroxyl Y39)
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Uniprot ID	P04908
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.