

Product Datasheet

Phospho-Smad2 (Ser250) monoclonal antibody (orb1723215)

Description	Phospho-Smad2 (Ser250) monoclonal antibody
Species/Host	Rabbit
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser250 of human Smad2
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 52 kDa; Observed MW: 58,62 kDa
Uniprot ID	Q15796
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.