

Product Datasheet

Myelin Basic Protein monoclonal antibody (orb1723854)

Description	Myelin Basic Protein monoclonal antibody
Species/Host	Rabbit
Reactivity	Hamster, Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human Myelin Basic Protein
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 33 kDa; Observed MW: 18 kDa
Uniprot ID	P02686
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.