

Product Datasheet

Phospho-LAT (Tyr200) monoclonal antibody (orb1723949)

Description	Phospho-LAT (Tyr200) monoclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Tyr200 of human LAT
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 28 kDa; Observed MW: 36-38 kDa
Uniprot ID	O43561
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.