

Product Datasheet

Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) polyclonal antibody (orb1725054)

Description	Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthesized peptide derived from human Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602)
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 64,80,75 kDa; Observed MW: 64,80,75 kDa
Uniprot ID	Q9P286, O96013, Q9NQU5
Dilution Range	WB: 1/500-1/1000 FC: 1/50-1/100
Expiration Date	12 months from date of receipt.