

Product Datasheet

MonoMethyl-Histone H3 (Lys18) polyclonal antibody (orb1725326)

Description	MonoMethyl-Histone H3 (Lys18) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	A synthesized peptide derived from human Histone H3 (mono methyl K18)
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Uniprot ID	P68431
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
Expiration Date	12 months from date of receipt.