

Product Datasheet

Phospho-alpha 1 Sodium Potassium ATPase (Ser23) polyclonal antibody (orb1725531)

Description

Phospho-alpha 1 Sodium

	polyclonal antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from rat ATP1 alpha1/Na+K+ ATPase1 around the phosphorylation site of Ser23. AA range:15-64
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 113 kDa; Observed MW: 113 kDa
Uniprot ID	P06685
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 ELISA: 1/10000
Expiration Date	12 months from date of receipt.