

## Product Datasheet

# Phospho-Connexin 43 (Ser368) polyclonal antibody (orb1725643)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | Phospho-Connexin 43 (Ser368)<br>polyclonal antibody                                      |
| <b>Species/Host</b>        | Rabbit                                                                                   |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | IHC-P, WB                                                                                |
| <b>Immunogen</b>           | Synthetic peptide of human GJA1                                                          |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.     |
| <b>Concentration</b>       | 1mg/ml                                                                                   |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                               |
| <b>MW</b>                  | Calculated MW: 43 kDa; Observed MW: 43 kDa                                               |
| <b>Uniprot ID</b>          | <b>P17302</b>                                                                            |
| <b>Dilution Range</b>      | WB: 1/500-1/1000 IHC: 1/50-1/100                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |