

Product Datasheet

14-3-3 polyclonal antibody (orb1725664)

Description	14-3-3 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 27-29 kDa; Observed MW: 27-29 kDa
Uniprot ID	P62258, Q04917, P31946, P27348, P31947, P61981, P63104
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.