

Product Datasheet

CD300b polyclonal antibody (orb1725784)

Description	CD300b polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Synthetic peptide from human protein at AA range: 61-110
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	A8K4G0
Dilution Range	IHC: 1/50-1/100 ELISA: 1/10000
Expiration Date	12 months from date of receipt.