

Product Datasheet

Histone H1 polyclonal antibody (orb1725829)

Description	Histone H1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Histone H1. AA range:1-50
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 22 kDa; Observed MW: 30 kDa
Uniprot ID	P16401, P16402, P10412
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 ELISA: 1/10000
Expiration Date	12 months from date of receipt.