

Product Datasheet

MYH polyclonal antibody (orb1725907)

Description	MYH polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MYH-pan around the non-acetylation site of Lys1394. AA range:1351-1400
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 223 kDa; Observed MW: 223 kDa
Uniprot ID	Q9UKX2 , Q9Y623 , P12882 , P11055 , P13533 , P13535 , P12883
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000
Expiration Date	12 months from date of receipt.