

Product Datasheet

CD1a (4C3) monoclonal antibody (orb1726030)

Description	CD1a (4C3) monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic peptide conjugated to KLH.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P06126
Dilution Range	IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.