

Product Datasheet

Phospho-Tyrosine (6G6) monoclonal antibody (orb1726033)

Description	Phospho-Tyrosine (6G6) monoclonal antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Purified Protein
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Refer to figures
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.