

Product Datasheet

beta III Tubulin (3H2) monoclonal antibody (orb1726060)

Description	beta III Tubulin (3H2) monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide of β III tubulin
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Uniprot ID	Q13509
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.