

Product Datasheet

MonoMethyl-Histone H3 (Lys9) (10F5) monoclonal antibody (orb1726077)

Description	MonoMethyl-Histone H3 (Lys9) (10F5) monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide of Histone H3 (Mono Methyl Lys9)
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Uniprot ID	P68431
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.