

Product Datasheet

MAP2K6 (8A11) monoclonal antibody (orb1726182)

Description	MAP2K6 (8A11) monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Purified recombinant human MKK6 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 37 kDa; Observed MW: 40 kDa
Uniprot ID	P52564
Dilution Range	WB: 1/500-1/1000 IP: 1/20
Expiration Date	12 months from date of receipt.