

Product Datasheet

SIRT6 (6C9) monoclonal antibody (orb1726259)

Description	SIRT6 (6C9) monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Purified recombinant human SIRT6 protein expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 39 kDa; Observed MW: 42 kDa
Uniprot ID	Q8N6T7
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Expiration Date	12 months from date of receipt.