

Product Datasheet

SENP2 Rabbit mAb Antibody (orb1726904)

Description	SENP2 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	Q9HC62
Expiration Date	12 months from date of receipt.