

Product Datasheet

GSTK1 Rabbit mAb Antibody (orb1727153)

Description	GSTK1 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	Q9Y2Q3
Expiration Date	12 months from date of receipt.