

Product Datasheet

Cyclin A1/A2 Rabbit mAb Antibody (orb1727556)

Description	Cyclin A1/A2 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P78396, P20248
Expiration Date	12 months from date of receipt.