

Product Datasheet

eIF4B (Phospho-Ser406) Rabbit mAb Antibody (orb1727675)

Description	eIF4B (Phospho-Ser406) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P23588
Expiration Date	12 months from date of receipt.