

Product Datasheet

MLKL (Phospho-Ser358) Rabbit mAb Antibody (orb1727702)

Description	MLKL (Phospho-Ser358) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	Q8NB16
Expiration Date	12 months from date of receipt.