

Product Datasheet

GFP Antibody (orb1727897)

Description	Chicken anti GFP [IQGFP6 Polyclonal]- Biotin
Species/Host	Gallus
Conjugation	Unconjugated
Immunogen	A synthetic peptide derived from the internal framework of green fluorescent protein
Target	GFP
Preservatives	Purified IgY supplied in Phosphate buffered saline (PBS) containing 0.01% Thiomersal and 1% trehalose, conjugated to Biotin
Concentration	1mg/ml
Note	For research use only
Isotype	IgY
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.