

Product Datasheet

Recombinant human TPOR protein, N-His (orb1727927)

Description	Recombinant human TPOR protein, N-His
Endotoxins	Not analyzed
Target	MPL
Preservatives	4M Urea
Form/Appearance	Lyophilized or Liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	N-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	297-482/635
Purity	> 90% as determined by SDS-PAGE
MW	27.4 kDa
Uniprot ID	P40238
Expiration Date	6 months from date of receipt.