

Product Datasheet

HLA-C Antibody (PE) (orb1727931)

Description	HLA-C Antibody (PE)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	Recombinant human HLA-C protein
Target	HLA-C
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	4 kDa
Uniprot ID	P10321
Expiration Date	12 months from date of receipt.