


Product Datasheet

Human MDGA1 ELISA Kit (orb1728063)

Description	Human MDGA1 ELISA Kit (96 Tests). Quantitate Human MDGA1 in cell culture supernatants, serum and pla...
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q19-G931
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q8NFP4
Sensitivity	30 pg/ml
Expiration Date	Please enquire.



Human
MDGA1
PicoKine
ELISA Kit
standard ...