


Product Datasheet

Human LRRTM3 ELISA Kit (orb1728067)

Description	Human LRRTM3 ELISA Kit (96 Tests). Quantitate Human LRRTM3 in cell culture supernatants, serum and p...
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E31-K419
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q86VH5
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Human
LRRTM3
PicoKine
ELISA Kit
standard...